


# NWF2046

## GYNECOLOGICAL TABLE



# NWF2046

## FEATURES AND HIGHLIGHTS



### Technical parameters

|                         |                 |
|-------------------------|-----------------|
| ▪ External size (LxWxH) | 1900x600x800 mm |
| ▪ Back plate up/down    | 7°/45°          |
| ▪ Seat plate up/down    | 8°/25°          |
| ▪ Leg plate up/down     | 25°/25°         |

### Technical configuration

|  |       |
|--|-------|
| <input checked="" type="checkbox"/> Crank Wheel              | 1set  |
| <input checked="" type="checkbox"/> Hand Bar                 | 1set  |
| <input checked="" type="checkbox"/> Anti-slip Bed Feet Cover | 4pcs  |
| <input checked="" type="checkbox"/> Mechanical Level         | 2sets |
| <input checked="" type="checkbox"/> Lithotomy Rod            | 1set  |
| <input checked="" type="checkbox"/> Shoulder Holder          | 1set  |
| <input checked="" type="checkbox"/> Stirrup                  | 1set  |
| <input checked="" type="checkbox"/> Arm Holder               | 1set  |
| <input checked="" type="checkbox"/> 3" Castors with Break    | 4pcs  |





### Structural design

The bed is composed of two parts: main bed and leg bed. The main bed table is divided into two parts: the back plate and the seat plate.



### Application

This delivery bed is suitable for general operations such as maternal delivery, prenatal examination, abortion surgery, gynecological examination, etc.




### Bed frame

The main bed and the auxiliary bed are separated separately, or they can be connected into one.



### Leg bed table

The upper and lower folds of the leg bed table are adjusted by the ratchet rack mechanism to ensure the comfort of the patient's legs during the operation and the main bed table folds to meet the needs of the posture of the operation.



### Mattress

Composed of 3 pieces, connected by a push button, with a sponge inside the artificial leather, leather fireproof, yellowing resistance, anti-mildew ball, twists and turns, abrasion resistance 100,000 times.



### Fixed and movable casters

4 non-slip foot pads. 4 movable casters, and The foot brake locking device is attached to ensure the stability of the leg bed when it is used, and it is more convenient to push.



### Anesthesia infusion stand

Use to hang medicine bottle during infusion, height can be adjusted.



### Shoulder Holder



### Arm Holder



### Lithotomy Rod



### Stirrup



### Hand Bar



# NUMBER WIN

MEDICAL TECHNOLOGY

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.



## 7x24

Online 24 Hours for Sales Service  
Timely Response for Clients' Problems

+86 18021231901

Suzhou Number Win Medical Technology Co.,Ltd.

No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

export@number-win.com

www.number-win.com



Company video



General catalogue